

LANTEX



SINGLE USE CIRCULAR STAPLER FOR RECTAL PROLAPSE AND HAEMORRHOIDS




SUZHOU LANTEX MEDICAL TECHNOLOGY COMPANY LIMITED

Building 7, 58 Xinting Road, High-tech Industrial Zone, 215151 Suzhou, China

Tel: +86 512 6616 7488

Fax: +86 512 6616 7488





PRODUCT SPECIFICATIONS

PRODUCT CODE	CARTRIDGE COLOR	STAPLE ROWS	OUTSIDE DIAMETER (mm)	CUTTING DIAMETER (mm)	STAPLE CROWN LENGTH (mm)	STAPLE LEG LENGTH (mm)	STAPLE COUNT (mm)	CLOSED STAPLE HEIGHT (mm)
CPH32	White	2	32.5	24	3.8	4.2	32	0.75-2.0
CPH34	Green	2	34.5	26	3.8	4.2	32	0.75-2.0
CPH36	Lime Green	2	36.5	28	3.8	4.2	36	0.75-2.0

Set contains: Circular Stapler for Prolapse and Haemorrhoids(CPH), Suture Threader, Anal Dilators (Round and "Butterfly"), Anal Dilator Inner Lining and Purse String Suture Anoscope