




PAJUNK®

*Modular
Laparoscopic Instruments
Flexible and safe*



Laparoscopic instruments

Variety in summary: Systematized modules

Always the right combination: As a modular system, the handles, the isolated shafts and

the instrument inserts may be combined in whichever size-orientated composition preferred.


Modules for PAJUNK® Laparoscopic-Instruments

Handles



Product	Item no.
Plastic handle for modular instrument inserts of the series 1292. Rotatable, with gold-plated HF connection. Automatic safeguard against distortion in the locked / engaged position.	1292-10-00
Plastic handle for modular instrument inserts of the series 1292. Rotatable, with gold-plated HF connection and ratchet. Automatic safeguard against distortion in the locked / engaged position.	1292-10-20
Plastic handle for modular instrument inserts of the series 1292. Rotatable, with gold-plated HF connection and switchable ratchet. Automatic safeguard against distortion in the locked / engaged position.	1292-10-30
Plastic handle for modular instrument inserts of the series 1292. Rotatable, with gold-plated HF connection and tilting ratchet underneath. Automatic safeguard against distortion in the locked / engaged position.	1292-10-40
Modular axial metal handle, with ratchet. Especially for use with needle holders 1292-95-xx. Without coagulation.	1292-20-30






















Insulated shafts




Product	Size	Working length	Item no.
Insulated shaft, adaptable to all handles and instrument inserts of the series 1292.	ø 5 mm	330 mm	1292-13-05
	ø 10 mm	330 mm	1292-13-10
	ø 5 mm	450 mm	1292-14-05
	ø 10 mm	450 mm	1292-14-10

Modular instrument inserts




Grasping forceps

	Product	Size	Working length	Item no.
	Dissection- and grasping forceps "DELPHIN", double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-11 1292-54-11
	Grasping forceps with reservoir, blunt, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-12 1292-54-12
	Grasping forceps with reservoir, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-13 1292-54-13
	Atraumatic grasping forceps, pointed, cross-cut serration, 18 mm, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-15 1292-54-15
	Mouse-tooth-type grasping forceps with 2 x 4 teeth, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-16 1292-54-16
	Alligator-type grasping forceps, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-17 1292-54-17
	Dissection- and grasping forceps, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-18 1292-54-18
	Dissection- and grasping forceps, fenestrated, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-19 1292-54-19
	Dissection- and grasping forceps "MARYLAND", double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-20 1292-54-20
	Dissection- and grasping forceps "MIXTER", double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-21 1292-54-21
	Grasping forceps "ALLIS", double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-22 1292-54-22
	Grasping forceps with fine cross-cut serration, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-23 1292-54-23
	Grasping forceps "BABCOCK", double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-26 1292-54-26
	Claw-type grasping forceps, 2 x 3 teeth, single action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-30 1292-54-30
	Grasping forceps "MAXI GRIP", double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-31 1292-54-31
	Grasping forceps "DORSEY", largely fenestrated, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-38 1292-54-38
	Atraumatic WAVE grasper, fenestrated, double action	ø 5 mm ø 5 mm	330 mm 450 mm	1292-53-39 1292-54-39
	Grasping forceps "MARYLAND", double action	ø 10 mm ø 10 mm	330 mm 450 mm	1292-63-20 1292-64-20
	Dissection- and grasping forceps "MIXTER", double action	ø 10 mm ø 10 mm	330 mm 450 mm	1292-63-21 1292-64-21
	Claw-type grasping forceps, 2 x 3 teeth, single action	ø 10 mm ø 10 mm	330 mm 450 mm	1292-63-30 1292-64-30
	Grasping forceps "BABCOCK", double action	ø 10 mm ø 10 mm	330 mm 450 mm	1292-63-31 1292-64-31




Biopsy Punch

	Product	Size	Working length	Item no.
	Biopsy punch with thorn, double action	ø 5 mm	330 mm	1292-73-14
		ø 5 mm	450 mm	1292-74-14




Clamps

	Product	Size	Working length	Item no.
	"BABCOCK" clamp, double action	ø 10 mm	330 mm	1292-83-10
		ø 10 mm	450 mm	1292-84-10
	"COLLIN" clamp, double action	ø 10 mm	330 mm	1292-83-15
		ø 10 mm	450 mm	1292-84-15
	Clamp with "DEBAKEY" teeth 2 x 3, jaw length 40 mm, double action	ø 10 mm	330 mm	1292-83-41
		ø 10 mm	450 mm	1292-84-41

Scissors

	Product	Size	Working length	Item no.
	Hook-type scissors, straight, single action	ø 5 mm	330 mm	1292-93-10
		ø 5 mm	450 mm	1292-94-10
	Dissecting scissors, straight, double action	ø 5 mm	330 mm	1292-93-13
		ø 5 mm	450 mm	1292-94-13
	Dissecting scissors, curved, "METZENBAUM", double action	ø 5 mm	330 mm	1292-93-14
		ø 5 mm	450 mm	1292-94-14












Needle holder

	Product	Size	Working length	Item no.
	Needle holder, straight, single action, with tungsten carbide coating and needle alignment	ø 5 mm	330 mm	1292-95-10
		ø 5 mm	450 mm	1292-96-10
	Needle holder, left curved, single action, with tungsten carbide insert	ø 5 mm	330 mm	1292-95-11
		ø 5 mm	450 mm	1292-96-11
	Needle holder, right curved, single action, with tungsten carbide insert	ø 5 mm	330 mm	1292-95-12
		ø 5 mm	450 mm	1292-96-12



Modular instrument shaft mountings

Series Ø 3 mm



Product	Working length	Item no.	Product	Working length	Item no.
 Dissection- and grasping forceps "DELPHIN", double action	250 mm 330 mm	1292-42-10 1292-43-10	 Dissecting scissors, straight, double action	250 mm 330 mm	1292-42-17 1292-43-17
 Atraumatic grasping forceps, pointed, cross-cut serration, double action	250 mm 330 mm	1292-42-13 1292-43-13	 Dissecting scissors, curved, "METZENBAUM", double action	250 mm 330 mm	1292-42-18 1292-43-18
 Universal grasping forceps, double action	250 mm 330 mm	1292-42-14 1292-43-14	 Hook-type scissors, straight, single action	250 mm 330 mm	1292-42-20 1292-43-20
 Dissection- and grasping forceps "MARYLAND", double action	250 mm 330 mm	1292-42-15 1292-43-15	 Micro-scissors, straight, single action	250 mm 330 mm	1292-42-21 1292-43-21
 Claw-type grasping forceps, 2 x 3 teeth, single action	250 mm 330 mm	1292-42-16 1292-43-16	 Biopsy punch, double action	250 mm 330 mm	1292-42-22 1292-43-22
			 Biopsy punch with thorn, double action	250 mm 330 mm	1292-42-23 1292-43-23

Adapters

Product	Item no.
 Adapter for high pressure rinsing, with female LL-connector for modular, nonseparable instrument mountings no. 1292-43-xx, as well as adaptable guidance tubes no. 1292-13-xx and 1292-14-xx	1292-43-98
 High pressure hose adapter, with male LL-connector	1292-43-99

PAJUNK GmbH
 Medizintechnologie
 Karl-Hall-Strasse 1
 D-78187 Geisingen / Germany
 Phone +49 (0) 77 04/92 91-0
 Fax +49 (0) 77 04/92 91-6 00
www.pajunk.com